


# Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)



Realy Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab) is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharynx.




25 Tests/kit

## Features:

- Specimen types: Nasal Aspirate Fluid/ Nasal Swab
- Testing time: 10 minutes
- Sensitivity: 96.17%
- Specificity: >99%

Product name	Ref.	Size	Storage Temp
Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)	K511416D	25 Test/kit	2-30°C
Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)	K511416D	5 Test/kit	2-30°C

## DIRECTIONS FOR USE:



HANGZHOU REALY TECH CO., LTD.  
4th Floor,#12 Building,Eastern Medicine Town,Xiasha  
Economic&Technology Development,310018 Hangzhou,  
Zhejiang,PEOPLE'S REPUBLIC OF CHINA  
<http://www.realytech.com>  
Email: info@realytech.com  
Tel: +86 0571-56050791/56050792/56050793

